

Small Universal Polarimeter

The Small Universal Polarimeter has circular scale of $\pm 180^\circ$ divided into full degrees. Together with the Vernier scale rotations can be measured to an accuracy of 0.1° and with a little practice they can be estimated to an accuracy of 0.05° .

The analyser is adjusted by rotating the Vernier scale, which, due to its large diameter, permits a most sensitive comparison of the half shadows. The scale, polariser, analyser and the channel for the samples are tilted 15° from horizontal, making operation more convenient and observation through the eye-piece more comfortable.

Polarising filters are used for the polariser and the analyser. The polariser is fitted with a thin two-piece quartz disc which produces a divided field of view, whose half shadow angle is dependent on the light used. With sodium light, $\lambda = 589 \text{ nm}$ it is 8° and with mercury light, $\lambda = 546 \text{ nm}$ it is 9.5° .

There are two possible light units. The less expensive halogen lamp for instructional purposes – it has a longer life and is ready for use immediately upon switching on – or the spectral lamp unit for higher absolute accuracy.

The lighting unit with halogen projector lamp is especially recommended for instructional purposes and for general measurements. Due to the low price of the unit and the long life of the halogen lamps, the unit is particularly economic. A further advantage is that the unit is ready for operation immediately upon switching on.

A transformer is needed to connect the unit to the normal mains supply.

Halogen lamps produce white light, and an interference filter is required to produce the monochromatic light of the desired wavelength. The filter is inserted into the light admission opening the sample channel. In this way measurements can be made at different wavelengths without time-consuming alterations to the apparatus (the halogen lamp does not need to be changed). Thus, for example, the nature of Optical Rotatory Dispersion can be very conveniently demonstrated.

Spectral lamps emit light of a defined wavelength (line spectrum). In order to block out undesired sidespectra, coloured glass filters should be used. Such unit may be used for any type of spectral lamp.

To connect spectral lamp lighting units to the mains, a current limiting ballast unit should be used. Such unit may be used for any type of spectral lamp.

Technical data:

Measurement range:	$\pm 180^\circ$
Resolution:	0.1°
Precision:	$\pm 0.1^\circ$
Light source:	halogen or spectral lamps
Wavelengths:	546 or 589 nm
Measuring tubes:	With bubble trap or filler cup, up to 200 mm
Mains adapter:	110/220 V / 50/60 Hz
Sensitivity:	up to 99,9% light absorption

SCHMIDT + HAENSCH GmbH & Co.

Waldstraße 80/81

D- 13403 Berlin

Germany

Phone: +49 30 / 41 70 72 – 0

Fax: +49 30 / 41 70 72 – 99

E-mail: sales@schmidt-haensch.de

www.schmidt-haensch.de

