

Washing water sterilizer 80 L

for Level 3/4/L3 safety labs



Hand washing, washing water sterilization, and disposal-all with one unit

Advantages

Hand washing basin with water supply and integrated sterilizer-all in one unit

Touch-free fittings provide well-tempered water

Sterilization starts and executes automatically

High degree of sterilization safety

Contaminated wash water cannot enter the waste water drain line

Only 600 mm wide-lab space is expensive

Washing water for 2 units operated in parallel: 60 litres per hour

Use

For washing and rinsing water in maximum containment genetic engineering labs and in other biological safety labs

Waste water sterilisation and disposal overnight or within 60 minutes

Steam Sterilization

Compared to chemical inactivation, steam sterilization offers:

- + maximum process safety
- + excellent reliability
- + environmentally friendly procedures
- + no health hazard to operators
- + low material and operating costs

VARIOKLAV®

●●● Washing Water Sterilizer