

Standard


CAL **Adjusting programme CAL**, for quick setting of the balance's accuracy. External test weights addit. price.

PCS **Parts counting:** reference quantities selectable. Display can be switched from piece to weight.

REZ **Net-total weighing:** weight of tare-cup and weight of components memorized in two separate stores.

IP 54 **Spray and dust protection:** protection class is indicated in the symbol.

230 V **Mains adapter:** 230 V / 50 Hz. Standard type D, on demand also available in type UK, AUS or USA.



DKD
Certificate
Possible

1
Warranty
Year

Compact Balance **CB**

For fast and precise weighing at a basic price



Large LCD display
with digit height 24 mm

Protective working cover
over user field and casing standard. Later order no. KERN CB-A 02

Battery supply possible
Block-battery (9 V) not standard. The balance has an AUTO-OFF function after 5 min to preserve the battery

Easy to use 4-key pad
with functions parts counting and mixing

Option

(Additional price)

AKU **Akku supply, external,** approx. 130 h. KERN KS-A 01.

DKD Calibration Certificates

DKD = German Calibration Service
KERN issues in its own accredited laboratories internationally valid DKD Calibration Certificates ...
■ for balances
■ for test weights

KERN is a manufacturer of test weights in all OIML classes of accuracy E₁, E₂, F₁, F₂, M₁-M₃, from 1 mg to 2000 kg

KERN - Your DKD Calibration Partner

Model	Read-out d	Weighing range Max	Repro-ducibility g	Linearity g	Min. piece weight at parts counting g	Net weight approx. kg
KERN						
CB 6 K 1	1	6	1	± 3	2	2.5
CB 12 K 2	2	12	2	± 6	4	2.5
CB 24 K 5	5	24	5	± 15	10	2.5

Total balance W x D x H 240 x 310 x 110 mm. Parts counting with reference 5, 10, 25, 50 pieces.

Your specialist dealer

