

Adjusting program (CAL): For quick setting of the balance's accuracy. External adjusting weight standard with KERN.



Data interface: RS 232 C, bi-directional, for connection of printer or PC.



Spray and dust protection IP 54: Protection class is indicated in the symbol to the left.



KERN PB 400-3 with cylindrical draught guard supplied as standard. Glas with metal top. Ø 150 mm, height 80 mm.



The special Tip:

New and innovative ceramic-quartz weighing system distinguished by high stability and precision.

Printer (additional price)

with data interface RS 232C for data transmission of weighing values. Paper width 55 mm. KERN 911-013.



KERN - Warranty: 3 years

DKD-Kalibrierlabor von KERN

(DKD = German Calibration Service)

Competence in weighing technology

- DKD-certificates for balances

Competence in adjust. weights

- KERN is a manufacturer of adjust. weights in all classes of accuracy E1 - M3, from 1 mg to 2 000 kg
- DKD-certificates for adjusting weights in the KERN laboratory



standard



Weighing range - Giant

Higher than standard weighing ranges are always advantageous, e.g. for heavy tare containers or simply as reserve. Nonetheless a reasonable price by focussing only the important functions.

Data-hold function: The display

is "frozen" via MODE key until the TARE key is pressed again.

Stand-by: No preheating required. Time saver!

Construction/design: Protective working cover supplied as standard. Compact metal casing, total ba-

lance WxDxH 170 x 270 x 75 mm. LED display, digit height 15 mm.

Common data: Stabilization time ≈ 2 sec, permissible ambient temperature +15 ... +30° C, mains adapter 230 V, net weight approx. 3.2 kg.

Model	Read-out	Weighing range	Reproducibility	Linearity	Weighing plate
	d	Max			stainl. steel
KERN	g	g	g	g	mm
PB 400-3	0.001	410	0.001	± 0.002	Ø 120
PB 4000-2	0.01	4100	0.01	± 0.02	Ø 160

Your specialist dealer:

Please request our balance catalogue!

KERN & Sohn GmbH, D-72322 Balingen, Postbox 40 52
 Tel. 00 49-[0]74 33-99 33-0 Fax 00 49-[0]74 33-99 33-149
 e-mail: info@kern-sohn.de Internet: www.kern-sohn.de

