


Standard

230 V **Mains adapter:** 230 V / 50 Hz. Standard type D, on demand also available in type UK, AUS or USA.

TAR **Taring:** place empty container on the pan, press TARE-key, zero digit appears. Repeatable → weigh in.

DKD
Certificate
possible



Standard

Warranty
Year

AKU
Option

Option

(Additional price)

AKU **Accu supply, external,** approx. 270 h. KERN KS-A 01.

Platform Balance **ECB**

- uncomplicated
- mobile
- low cost



Weighing plate
stainless steel
320 x 260 mm

Extremely flat construction

Only 60 mm height

Large LCD display
with digit height 14 mm

Fast

Stabilisation time
only approx. 2 sec

Easy to use

two keys only



This balance is especially suitable for mobile use thanks to its battery supply and its flat construction



Battery supply possible

Batteries size AA (6 x 1.5 V) not standard. The balance has an AUTO-OFF function after 3 min to preserve the battery. Service life approx. 100 h

Model	Read-out	Weighing range	Reproducibility	Linearity	Net weight approx.
KERN	d g	Max kg	g	g	kg
ECB 20 K 20	20	20	20	± 40	1.7
ECB 50 K 50	50	50	50	± 100	1.7

Total balance W x D x H 320 x 300 x 60 mm.

DKD Calibration Certificates

DKD = German Calibration Service

KERN issues in its own accredited laboratories internationally valid DKD Calibration Certificates ...

- for balances
- for test weights

KERN is a manufacturer of test weights in all OIML classes of accuracy E₁, E₂, F₁, F₂, M₁-M₃, from 1 mg to 2000 kg

KERN - Your DKD Calibration Partner

Your specialist dealer

KERN & Sohn GmbH, D-72322 Balingen, Postbox 4052
Tel. 00 49-(0)74 33-99 33-0 Fax 00 49-(0)74 33-99 33-149
e-mail: info@kern-sohn.de Internet: www.kern-sohn.com

