

Standard

TAR **Taring:** place empty container on the pan, press TARE-key, zero digit appears. Repeatable → weigh in.

TAR
Standard

230 V
Option


1
Warranty
Year

Option

(Additional price)

230 V **Mains adapter 230 V / 50 Hz, AUTO-OFF** remains active, KERN 440-902.

Battery Model EMB
Budget price, with large weighing plate



Socket
standard for mains adapter

CE

Weigh in function
simplifies mixing

Large LCD display
with digit height 16 mm

Model	Read-out d	Weighing range Max g
KERN	g	g
EMB 200-1	0.1	200
EMB 2000-0	1	2000
EMB 5 K 5	5	5000

Total balance W x D x H 170 x 240 x 38 mm.
Net weight approx. 600 g.

Ready for use

Batteries size AA (4 x 1.5 V) inclusive. The balance has an AUTO-OFF function after 3 min to preserve the battery



DKD Calibration Certificates

DKD = German Calibration Service

KERN issues in its own accredited laboratories internationally valid DKD Calibration Certificates ...

- for balances
- for test weights

KERN is a manufacturer of test weights in all OIML classes of accuracy E₁, E₂, F₁, F₂, M₁-M₃, from 1 mg to 2000 kg

KERN - Your DKD Calibration Partner

Your specialist dealer

KERN & Sohn GmbH, D-72322 Balingen, Postbox 4052
Tel. 00 49-(0)74 33-99 33-0 Fax 00 49-(0)74 33-99 33-149
e-mail: info@kern-sohn.de Internet: www.kern-sohn.com

