



MSZ 5000

Stereo Zoom Microscope-Serie MSZ 5000

Zoom - the adjustable alternative, with accessories up to 180x magnification

[Info](#) [Order](#)

(Picture top)

Stereo Zoom Microscope without illumination.

- ▶ 10x plano eyepieces
- ▶ with diopter compensation
- ▶ adjustable interocular distance
- ▶ Zoom lense from 0.7x to 4.0x, for zoom magnification from 7x to 40x (extended up to 160x)
- ▶ Metal stand with working distance 70 mm



Accessories:

- ▶ Eyepieces 20x (14x - 80x magnification)
- ▶ Additional lens 0.5x (3,5x - 20x magnification)
- ▶ Additional lens 2x (14x - 80x mag) with 20x eyepieces -160x



MSZKL 5000

- ▶ The MSZKL 5000 is equipped as MSZ 5000, but with built-in cold light source. Ideal for incident and transmitted light microscopy.

MSZ 5500 Stereo Zoom Microscope

- ▶ Like MSZ 5000, but with an additional photo tube for video or film camera adaption.

MSZ 5400 Stereo Zoom Microscope (Picture Bottom)

- ▶ Stereo Zoom Microscope with built-in illumination.



- ▶ 10x plano eyepieces, with diopter compensation and adjustable interocular distance
- ▶ Zoom lenses from 0.7x to 4.5x, for zoom magnification from 7x to 45x (extended up to 180x)
- ▶ Metal stand with 90 mm working distance
- ▶ Built-in illumination for direct and transmitted light (6 V, 15 W)

Accessory:

- ▶ Eyepieces 20x 14x - 90x magnification
- ▶ Additional lens 0.5x 3.5x - 22.5x magnification
- ▶ Additional lens 2x 24x - 180x magnification

**5600 Stereo Zoom Microscope**

- ▶ Same as MSZ 5400, but supplied with additional photo tube for video or film camera adaption.

[Addresses & Contact](#)[Information/Offer Order](#)[Articles & Papers](#)[History](#)[Home](#)

A. Krüss Optronik GmbH · Alsterdorfer Strasse 220 · D-22297 Hamburg

Tel.: +49 40 514 31 70 · Fax: +49 40 51 25 22 · E-Mails: info@kruess.com / [Webmaster](#)

Rua Pereira de Almeida, 31 Pça da Bandeira
Rio de Janeiro – RJ. CEP. 20260-100
Tel. (21) 502-0108; 502-2978; 273-3667
Fax.(21) 504-0868
E-mail info@precitech.net
HP: www.precitech.net